


Transmission	T(avg) > 90% (505.5 - 544.5 nm)
Blocking	OD(avg) > 5 (275 - 489 nm) OD(avg) > 3.3 (489 - 494 nm) OD(avg) > 3.3 (560 - 565 nm) OD(avg) > 5 (565 - 740 nm)
Polarization	Average
Angle of Incidence	0.0° (± 5.0°)
Cone Half Angle	0.0° (± 0.0°)
AR Coating	None
Filter Housing	Mounted in a black anodized aluminum housing
Housing Thickness	3.50 mm (± 0.10 mm)
Filter Diameter	25.00 mm (+ 0.00 / - 0.10 mm)
Clear Aperture (CA)	≥ 22 mm
Substrate Thickness	2.00 mm (± 0.10 mm)
Wedge	≤ 20 arcseconds
Flatness	No Specification
Substrate Material	Borofloat
Transmitted Wavefront Error	≤ 1 wave PV per inch @ 633 nm
Scratch / Dig	60/40
Bevel	≤ 0.3 mm @ 45° (as necessary)
Edge Chip	As per ANSI/OEOSC OP1.002-2009, American Standard
Filter Marking Text	BrightLine Basic™ Fluorescence Filter 525/39
Packaging	Filters will be individually packaged.
Package Labeling	Packaging will be labeled with Semrock logo, part number, revision, lot code, and work order number.
Special Requirements	100

Drawing below not to scale - for Reference Only



- Optical density specifications greater than 5.5 OD guaranteed by design.
- RoHS compliant
- Designed to meet or exceed the requirements for environmental and physical durability set forth in U.S. Military specifications MIL-STD-810F, MIL-C-48497A, MIL-C-675C, as well as the international standard ISO 9022-2

	Drawn By: A. Hill	Date: July 27, 2010	Scale: None	Page 1 of 1	
This document contains proprietary information of Semrock, Inc. and is to be retained in confidence not to be reproduced or copied in whole or in part and not to be used or incorporated in any manner, except under expressed written agreement with Semrock, Inc.	Semrock, Inc. - A Unit of IDEX Corp. 3625 Buffalo Road, Suite 6 Rochester, NY 14624 USA		Part / Dwg Number		Rev
			FF01-525/39-25		INIT
			525/39 BrightLine® Bandpass Filter, 25 mm		
	+1-585-594-7000 www.semrock.com semrock@idexcorp.com				